INSTRUCTIONS FOR USE OF:

IM3[®] EVOLUTION RANGE

THE VETERINARY DENTAL UNITS

ENGLISH MANUAL

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ČESKÝ MANUÁL

Models **Evolve/Advanced/Professional**

Version 1.14092023

www.im3vet.com

CONTENTS

1. INTRODUCTION	4
» About this document	4
» Definition of symbols	5
› WARNINGS SYMBOLS	5
» Definition of symbols	6
› Other symbols used	6
» Responsibility of the manufacturer	7
» Explanation of warranty terms	7
» Copyright information	7
2. INSTALLATION	8
» Evolution Cart with Compressor	8
» Operating & Transport conditions	9
» Install requirements	9
» Unboxing the unit	10
3. OPERATION - UNIT FEATURES	11
» SWITCH the unit ON	15
» Instrument bridge in Start Position	15
» Foot Control	16
» Instrument control	16
» 3-1 Syringe	16
» Display Screen	17
› Intended Use	17
› Overview of the Display	17
› Factory Settings - RESTORATION & PROSTHETICS	18
› Factory Settings - PROPHYLAI	18
› Operation - CHANGING AND SAVING SETTINGS	19
› Operation - FORWARD / REVERSE SETTINGS	19
› Operation - ENDODONTICS SETTINGS	20
› Operation - GENERAL SET-UP MENU	20
› Standby Mode	20
› Operation - BACK TO ORIGINAL PROGRAM	20
› Operation - RESET FACTORY SETTINGS	20
› RESTORATION & PROSTHETICS	21
› ERROR MESSAGES	22
› Technical Data	23
» Electric Micromotor iM3 em-12 I ga	24
» Air Instrument	24
› Low & High Speed Handpiece Pressure	24
› Ultrasonic scaler	25
› Overview of Display	26
› ERROR Messages for Scaler if using the display	26
› Changing / Saving the power on the Ultrasonic Scaler	26
› Suction	27
› Distilled Water Bottle & CLS Bottle	28
› Bur Replacement - iPush Tool	29
› iM3 Advantage Low Speed Handpiece	29
› iM3 High Speed Handpiece	30
» iM3 Bottle System	31
› Installation & Operating Instructions	31
› 1. Safety	32
› 2. Accessories	34
› 3. Technical data	35
› 4. Function	36
› 5. Handling	37
› 6. Disinfection and cleaning	37
› 7. Reprocessing	39
› 8. Maintenance	42
› 9. Tips for operators and service technicians	43
› 10. Manufacturer of Bottles	43
4. INFECTION CONTROL	44
» cross-contamination	44
› cleaning, disinfection & sterilization procedures	44
» General Disinfection of unit surfaces, handles and skirt	45
» Infection control of the Silicon mat	45
» Disinfection of instrument hoses	45
» Infection control of other instruments	46
» Disinfection of unit water lines	46
» High speed handpiece	47
» Low speed handpiece	47
5. UNIT MAINTENANCE	48
› Technical Data Unit	48
› Display	48
› Ultrasonic Scaler Handpiece	48
› Micromotor	49
› Compressor	49
6. UNIT ANNUAL SERVICE	50
› Evolve Dental Unit	50
7. LEGAL INFORMATION	51
› EMC - Electromagnetic compatibility	51
› CE - European conformity	51
› UKCA - UK Conformity Assessed	51
8. SPARE PART & SERVICE LIST	52
» AS001- Advanced service kit used Evolve	52
» AS002- Advanced service suction kit used on Advanced and Professional	52

CONTENTS cont'd

9. TROUBLE SHOOTING & MAINTENANCE 53

- » **Compressor 53**
 - › Identifying leaks 53
 - › Adjusting Regulator Air Pressure 54
 - › Compressor Maintenance 54
 - › Handpiece Maintenance 55
 - › Handpiece Trouble Shooting 56
 - › Low Speed Trouble Shooting 57
 - › Cartridge Toggle 57
 - › Ultrasonic Scaler Handpiece 57

10. TECHNICAL PRODUCT DETAILS 58

1. INTRODUCTION

ABOUT THIS DOCUMENT

Please read this manual carefully and explore all the equipment's extraordinary features:

- ✓ The unit shall be installed as described in **Section 2**.
- ✓ **Section 3** describes how to use the equipment.
- ✓ Read details about infection control and cleaning in **Section 4**.
- ✓ Maintain the unit as described in **Section 5**.
- ✓ The Evolution unit must be inspected and serviced every 12 months by an iM3 authorized service provider to ensure safe operation– see details in **Section 6**.
- ✓ In **Section 7** we list important legal information.
- ✓ Please see **Section 8** for a list of symbols used and **Section 11** provides technical product details.
- ✓ **Section 9** contains a complete list of accessories, detachable parts and consumables.
- ✓ **Section 10** contains the Evolution trouble shooting and daily weekly maintenance guide.

For more information, please contact customer service in your country :

EU & UK: support@im3dental.com | US: info@im3usa.com | AU: sales@im3vet.com

or you can contact your local distributor.

These installation and operating instructions form part of the unit. Failure to comply with the instructions and information in these installation and operating instructions means that iM3 will not be able to offer any warranty or assume any liability for the safe operation and the safe functioning of the unit.

WARNING! No unauthorized installation, service or modification of this equipment is allowed. Unauthorised modification could result in a malfunction or physical injury!

DEFINITION OF SYMBOLS

Please see below a list of symbols which can be found on:

- » instructions for use manual
- » the device/unit
- » the packaging

WARNINGS SYMBOLS

The warnings in this document are intended to draw your attention to possible injury to persons or damage to machinery.

The warnings are structured as follows:

 **SIGNAL WORD:** Description of the type and source of danger
 The possible consequences of ignoring the warning are described here.
 Follow these measures to avoid the danger.

The signal word differentiates between four levels of danger:

- » **DANGER**
Immediate danger of severe injury or death
- » **WARNING**
Possible danger of severe injury or death
- » **CAUTION**
Risk of minor injuries
- » **NOTICE**
Risk of extensive material/property damage

The following warning symbols are used:

	WARNING: General Warning Symbol
	WARNING: Risk of injury
	ATTENTION : something could be damaged.
	PROTECTION : Bio-Hazard warning
	PROTECTION : Wear protective gloves
	PROTECTION : Wear protective goggles
	PROTECTION : Use mouth & nose protection
	PROTECTION : Wear protective clothing

 These installation and operating instructions form part of the unit. Failure to comply with the instructions and information in these installation and operating instructions means that iM3 will not be able to offer any warranty or assume any liability for the safe operation and the safe functioning of the unit.



DEFINITION OF SYMBOLS

Please see below a list of symbols which can be found on:

- » instructions for use manual
- » the device/unit
- » the packaging

OTHER SYMBOLS USED

	General information		Important information
	Consult Instructions for Use		Manufacturer
	Date of Manufacture		Expiry Date
	Product Code		Serial Number
	Medical Device		Data structure in accordance with Health Industry Bar Code
	CE labelling: compliance with EU legislation		UKCA: UK Conformity Assessed for products sold within Great Britain.
	UL Component Recognition Mark indicates compliance with Canadian and U.S. requirements		Data Matrix Code for product information including Unique Device Identification
	Caution! According to Federal law, this medical device may only be sold by or on the order of a dentist, physician or any other medical practitioner licensed by the law of the state in which he or she practices and intends to use or order the use of this medical device.		
	Sterilizable up to the stated temperature		Thermo washer disinfectable
	Humidity limitation		Electrostatic sensitive devices
	Permitted temperature range		Supply voltage of the device
	Type B applied part (not suitable for intracardiac application)		Suitable for people with pacemakers or implanted defibrillators
	Fragile, handle with care		This way up
	DO NOT dispose of with domestic waste		Thread system: Q-Link Ultra style connection P6 style connection

These installation and operating instructions form part of the unit. Failure to comply with the instructions and information in these installation and operating instructions means that IM3 will not be able to offer any warranty or assume any liability for the safe operation and the safe functioning of the unit.

RESPONSIBILITY OF THE MANUFACTURER

The manufacturer can only accept responsibility for the performance of the device when it is used in compliance with the following directions:

- » The device must be used in accordance with these **Instructions for Use**.
- » Modifications or repairs must **ONLY** be undertaken by iM3 or an authorised iM3 service partner.
- » The device has no components that can be repaired by the user.
- » The electrical installation at the premises must comply with the regulations laid out in **IEC 60364-7-710 (Installation of electrical equipment in rooms used for approved purposes)** or with the regulations applicable in your country.
- » Unauthorised opening of the device invalidates all claims under warranty and any other claims.
- » Improper use, unauthorised assembly, modification or repair to the device, non-compliance with our instructions or the use of accessories and spare parts which are not approved by iM3, invalidates all claims under warranty and any other claims.
- » The appliance is designed to be used by a qualified Veterinary professional **ONLY**. The appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- » Children should be supervised to ensure that they do **NOT** play with the appliance.
- » **FOR VETERINARY DENTAL USE ONLY.**

EXPLANATION OF WARRANTY TERMS

- » **This iM3 device has been manufactured with great care by a qualified team member. A wide variety of tests and controls ensures faultless operation. Please note claims under warranty can only be validated when all the directions in the instructions for use have been followed.**
- » As the manufacturer, iM3 is liable for material or manufacturing defects within the 24 months from the date of purchase. Accessories and consumables (tips, tip changer, nozzle cleaner, pump cartridge, coolant hose, coolant filter, straw, O Rings) are excluded from the warranty.
- » We accept no responsibility for damage caused by incorrect handling or by repairs carried out by third party not authorised to do so by iM3.
- » Claims under warranty accompanied by proof of purchase, must be sent to the vendor or to an authorised iM3 service partner.
- » The provision of service under warranty extends neither the warranty period nor any guarantee period, repair and returns.
- » In the event of operating malfunctions immediately contact an authorised iM3 service partner.
- » Repairs and maintenance work must only be undertaken by an authorised iM3 service partner.

For a list of iM3 Authorised Service Partners in your country please visit:

Europe & UK: www.im3vet.eu

USA: www.im3vet.com

Australia/Asia: www.im3vet.com.au

COPYRIGHT INFORMATION

All circuits, processes, names, software programs and units mentioned in this document are protected by copyright. The installation and operating instructions may not be copied or reprinted, either in full or in part, without written authorisation from iM3.

2. INSTALLATION

EVOLUTION CART WITH COMPRESSOR


Fig. 1 - Top Tray




Fig. 2 - Bottom of Unit

Fig. 3 - Full height of Evolve Dental Unit

Wall Mounted Version

Fig. 4 - Wall connection

OPERATING & TRANSPORT CONDITIONS

CONDITION	OPERATION	TRANSPORT & STORAGE
Temperature	+10°C to + 35°C	- 20°C to + 70°C
Relative humidity	15% to 80%	8% to 80%
Air Pressure	800 hPa to 1060 hPa	700 hPa to 1060 hPa
Installation Altitude	Max 2,000 metres above sea level	

INSTALL REQUIREMENTS

ROOM	REQUIREMENT WITH TABLE
Width	200cm
Length	300cm
Height	> 200cm

ELECTRICAL	REQUIREMENTS
Mains Voltage	110V or 240V (at Option) Non switchable
Mains Frequency	50/60Hz
Primary Fuse	10A(110V) 5A (240V) Circuit Breaker

WATER	REQUIREMENTS
Type	Distilled or Ionized water only to be used
Filter	iM3 Straw fitted to water bottle, (replaced yearly)
Extras	iM3 CLS present in CLS bottle or iM3 straw fitted, (replaced yearly)

The Evolution Dental unit must be installed in accordance with the instructions/unboxing video or by an iM3 Authorized Service Provider.
Wall mounted version needs to be installed only by an iM3 Authorised Service Provider.

WARNING: To avoid risk of electrical shock, this equipment must be connected to a supply with protective earth.

To avoid the risk of electric shock always switch off the power to the unit in the electrical panel before opening or touching internal components.


Unboxing instructions for the Evolution range can be downloaded from www.im3vet.com


UNBOXING THE UNIT

1


Cut the straps holding the box to the pallet.

2


Remove the lid from the box and the top cardboard spacer.

3


Remove individual boxes in numerical order, by lifting upwards using provided finger holes.

4


Keep the numbered boxes as they contain handpieces and other items needed for your dental unit.

5


With the boxes removed, lift up the outer box, being careful not to hit the display screen if installed.

6


Remove inner box which surrounds your new Dental Unit.

7

Pull down the box side lip, unlock the wheels, hold the metal located under the blue skirt at the base of the unit and roll out - as above.
DON'T use blue handles or corian top to roll out

8

Install all other additional elements as per instructions supplied within the box. See note below

Fig.5 - Step-by-Step unboxing instructions

DO NOT lift the unit by any plastic part of the unit or by the Corian Top. Only lift from the metal sub structure under the skirt.

- » Please install handpieces as per instructions placed within the box.
- » When installing the straw, leave for 24hrs to clean internally before using the unit for the first time. Follow the straw instructions fully before use.
- » Fill the water bottles with distilled water. Fill the CLS bottle with distilled water and CLS as per the mixing instructions printed on the CLS bottle. Switch the unit on and you are now ready to use the machine once the instructions in section 3 Operation, have been read.

3. OPERATION - UNIT FEATURES


Fig.6a - UNITS FRONT FEATURES

i A complete technical description of the Evolution Range is available at: www.im3vet.com

! To avoid injury to persons or material, do not use the Evolution Range or its accessories if signs of operational, electrical or mechanical defects are found.

i Please note that the Evolution Range must be operated in accordance with the manual and by Veterinary Professionals Only.

! **!** Exercise caution when using the unit in combination with other equipment that can move.

! **!** Use of non iM3 Equipment adjacent to, or placed on, this equipment should be avoided because it could result in improper use.


Fig.6a -UNITS INSTRUMENT BRIDGE

A complete technical description of the Evolution Range is available at: www.im3vet.com

To avoid injury to persons or material, do not use the Evolution Range or its accessories if signs of operational, electrical or mechanical defects are found.

Please note that the Evolution Range must be operated in accordance with the manual and by Veterinary Professionals Only.

Exercise caution when using the unit in combination with other equipment that can move.

Use of non iM3 Equipment adjacent to, or placed on, this equipment should be avoided because it could result in improper use.


Fig.6c -UNITS SIDE LOCATED FEATURES

A complete technical description of the Evolution Range is available at: www.im3vet.com

To avoid injury to persons or material, do not use the Evolution Range or its accessories if signs of operational, electrical or mechanical defects are found.

Please note that the Evolution Range must be operated in accordance with the manual and by Veterinary Professionals Only.

Exercise caution when using the unit in combination with other equipment that can move.

Use of non iM3 Equipment adjacent to, or placed on, this equipment should be avoided because it could result in improper use.


Fig.6d - UNITS BACK LOCATED FEATURES

A complete technical description of the Evolution Range is available at: www.im3vet.com

To avoid injury to persons or material, do not use the Evolution Range or its accessories if signs of operational, electrical or mechanical defects are found.

Please note that the Evolution Range must be operated in accordance with the manual and by Veterinary Professionals Only.

Exercise caution when using the unit in combination with other equipment that can move.

Use of non iM3 Equipment adjacent to, or placed on, this equipment should be avoided because it could result in improper use.

SWITCH THE UNIT ON

Switch the unit ON (and OFF) using the main switch (see Fig. 7)

The Unit is ready for use after the compressor has reached pressure or the monitor has switched ON, whichever is last.

The compressor motor will switch OFF when the air tank is filled. The compressor motor will switch ON and OFF automatically as more air is needed.

Incase of emergency, use this switch to turn Off the unit.

Fig. 7 - Switch (ON/OFF) located at side of unit

INSTRUMENT BRIDGE IN START POSITION

Always place the handpieces in the start position in the instrument bridge when placing the animal on and off the table as shown below.

All instruments should be placed into the instrument holders in the instrument bridge (as pictured Fig. 8) when not in use. This reduces the likelihood of the instruments being damaged if they are knocked off the table and hit the ground.

Instrument Bridge from left to right (fig.8)

- a. Suction
- b. Ultrasonic Scaler
- c. 3-1 air water syringe
- d. Low Speed Air (air instrument)
- e. High Speed Air (air instrument)
- f. Electric Micromotor

Please refer to the supplied instrument manuals that comes with each handpiece for relevant information on your instruments

Fig. 8 - Instrument Bridge
(L to R: Suction, Ultrasonic Scaler, 3-1 air Water Syringe, Low Speed Air, High Speed Air, Electric Micromotor)

FOOT CONTROL

The Foot Control operates the handpieces (see Fig.9). Pressing down onto the foot control will start the selected instrument. Lifting off the foot control will stop the instrument.

The foot control should be stored on the dental unit shelf when not in use (see Fig. 10 below).

Fig. 9 - Evolve's Foot Control


Fig.10 - Storage of Foot control on Unit

INSTRUMENT CONTROL

The position of each instrument is listed in *figure 8* on the previous page, and this is classed as the *Start Position* of your instruments.

Activate the first instrument by lifting the handpiece up and pressing the foot control.

To avoid unintentional activation of the second instrument the first instrument needs to be placed back into the holder in the start position before the second instrument can be used.

The 3-1 syringe and suction can be used in conjunction with any other instrument.

Fig.11 - instruments in Start Position

3-1 SYRINGE

The 3-1 Syringe is operated by pressing the buttons. When pressed, the left button provides water, and the right button provides air. Pressing both buttons simultaneously will deliver a mist.

Fig.12 - 3-1 Syringe

DISPLAY SCREEN

INTENDED USE

Component for displaying and setting parameter of controlling devices.

Fig.13 Units Display Screen

OVERVIEW OF THE DISPLAY


Fig.14 - Display Screen's Operational overview

FACTORY SETTINGS - RESTORATION & PROSTHETICS

- P** Program 1 – Display shows P1: forward drive, 40,000 rpm
- P** Program 2 – Display shows P2: forward drive, 20,000 rpm
- P** Program 3 – Display shows P3: forward drive, 2,000 rpm

Note: the transmission / reduction ratio of the transmission instrument, speed of the instrument (e.g. rotary instrument) depends on the handpiece used.

The Display shows the speed of the motor

FACTORY SETTINGS - PROPHYLAXIS

- P** Program 1, display shows P1: 30
- P** Program 2, display shows P2: 20
- P** Program 3, display shows P3: 15

FACTORY SETTINGS - ENDODONTICS

(identical for 2:1 and 16:1 reduction ratio)

- P** Program 1, display shows P1: 1:1, 300 rpm, 3.0 Ncm
- P** Program 2, display shows P2: 1:1, 300 rpm, 2.5 Ncm
- P** Program 3, display shows P3: 1:1, 300 rpm, 1.5 Ncm
- P** Program 4, display shows P4: 1:1, 300 rpm, 1.0 Ncm
- P** Program 5, display shows P5: 1:1, 300 rpm, 0.5 Ncm

The display shows the speed of the root canal instrument

FACTORY SETTINGS - KEYS

PLUS key:
Increase speed / power, possible during use

MINUS key:
Reduce speed / power, possible during use
By holding PLUS / MINUS key the values are continuously increased / decreased

PLUS Ncm key:
Increase torque

MINUS Ncm key:
Reduce torque

ENDO key:
Switch between Restoration & Prosthetics operation, and Endodontics operation

FORWARD / REVERSE operation key:
Changing from forward to reverse operation is not possible during use.
[inactive in Endodontics / Prophylaxis operation]

OPERATION - CHANGING AND SAVING SETTINGS

i The speed / power can be set at P for minimum 2,000 rpm to maximum 40,000 rpm / 1 to maximum 40

P Step 1: Press the program key [P]

+ Step 3: Increase Speed

+ Step 2: Increase Torque
[inactive in Restoration & Prosthetics operation]

- Step 4: Reduce Speed

+ Step 3: Reduce Torque
[inactive in Restoration & Prosthetics operation]

P Step 5: Press the program key [P] for approximately 2 seconds to save.

i The set values flash and an acoustic signal is emitted to confirm that all the settings have been saved.

OPERATION - FORWARD / REVERSE SETTINGS

i Factory Setting – Forward Operation.

↻ Step 1: Press the forward / reverse operation key

i »!« and ↻ appear on the display and an acoustic signal is emitted.

i »!« and ↻ will flash continuously. An acoustic signal is emitted three times in succession before the motor is started in reverse operation, via the foot control.

P Step 2: Press the program key [P] for approximately 2 seconds to save.

i If the program settings are not saved in reverse operation, the device will switch back to Forward Operation automatically, when you change programs.

OPERATION - ENDODONTICS SETTINGS

Step 1: Press the Endo key to switch between Restoration & Prosthetics operation, and Endodontics operation.

Endo key Flashes: »!«, » **turn off**«,  appear on the display.

Step 2: Switch off Coolant supply at the dental unit.

Step 3: Press program key [P] to select programs [1 – 5].

STANDBY MODE

The Standby Mode is automatically activated after 3 minutes. You can exit standby mode by actuating the foot control or pressing the keys.

OPERATION - GENERAL SET-UP MENU

You can return to the set-up menu from any program by pressing on the PLUS and MINUS keys at the same time.

Step 1: Press PLUS and MINUS keys simultaneously for approximately 3 seconds. »Setup« appears on the display.

Step 2: Press program key P to navigate in the set-up menu.

Step 3: Press the PLUS / MINUS key to change the settings in the respective set-up menu item.

Step 4: Press program key P for approximately 2 seconds to save.

The set values flash and an acoustic signal is emitted to confirm that all the settings have been saved.

OPERATION - BACK TO ORIGINAL PROGRAM


Press PLUS and MINUS keys simultaneously for approximately 3 seconds to exit the Set-up Menu.

The original program appears on the display.

OPERATION - RESET FACTORY SETTINGS

Simultaneously press the program P, Endo and – Ncm keys, for approximately 3 seconds, to return the current program to factory defaults.

RESTORATION & PROSTHETICS

»Pedal«

- > ON / OFF
- > STEP: variable from 2,000 rpm to value set [maximum 40,000 rpm]
Endodontics: variable = maximum speed
Prophylaxis: power variable from 1 to maximum 40

»Sound«

- > ON
- > OFF

»LED«

 >Setting for fade-out time: from 0 to maximum 60 seconds.
Factory default: 5 programs

»P1 – P5«

 >Set the number of programs to be available in Restoration and Prosthetics operation, Factory default: 3 programs

»Speed«

- >rpm =Absolute speed in rpm
- >% = speed displayed as a percentage

P1	100%
P2	50%
P3	5%

Speed rpm	%
40,000	100
30,000	75
20,000	50
10,000	25
4,000	10
2,000	5

»1:1 / 2:1«

	Speed rpm		Torque Ncm	
	Max.	Min.	Max.	Min.
WH-56 1:1	3,000	100	3	0.2
EB-79 [2:1]	3,000	100	5	0.2

»AR – AF«

- > Auto Forward / Auto Reverse:
When the preset torque is reached, the handpiece drive runs in an anti-clockwise direction. The drive then stops, freeing the file.

»P1 – P9«

- > Saving additional programs / root canal instruments

ERROR MESSAGES

DISPLAY	DESCRIPTION OF ERROR	SOLUTIONS
Error 1	Overheating / overloading of electronics	> ❶ Disconnect device from power supply > ❷ Wait 5 minutes and allow the system to cool > ❸ Switch the device back on again and restart the function
Error 2	Pedal of the foot control pressed during switching on	> Do not press the pedal of the foot control
Error 3		
Error 4	Display keys pressed during switching on	> Do not press the display keys
Error 5	Running time limiter as a result of 15 minutes of continuous operation	> Check the pedal of the foot control > Do not press the pedal of the foot control any more [release completely]
Error 6	Error with »Electric motor / Endo contra-angle handpiece / Piezo scaler«	> ❶ Check the electric motor / Endo contra-angle handpiece / scaler is correctly attached to the coupling > ❷ electric motor / Endo contra-angle handpiece / replace tip
Reboot		> Switch the power supply off and on again > Restart system
Error	e.g.: 05 6303	> Contact an authorised iM3 service partner immediately

If the error messages described cannot be resolved a check by an authorised service partner is required.

For a list of iM3 Authorised Service Partners in your country please visit:

Europe & UK: www.im3vet.com

USA: www.im3vet.com

Australia/Asia: www.im3vet.com.au

TECHNICAL DATA

MODEL	D-2.0 E
Power supply:	5 V DC $\pm 10\%$, from W&H controller $\sqrt{\text{---}}$ (supply voltage of the device)
Power consumption:	1.5 VA
Dimensions in mm [WxDxH]:	155.5 x 29 x 122
Weight:	280 g

✓ AMBIENT CONDITIONS

Temperature during storage and transport	-40°C to +70°C [-40°F to +158°F]	
Humidity during storage and transport	8% to 80% [relative], non-condensing	
Temperature during operation	+10°C to +35°C [+50°F to +95°F]	
Humidity during operation	15% to 80% [relative], non-condensing	
Permitted ambient pressure	70 – 106 kPa	
Pollution level	2	
Overvoltage category	II	
Altitude	up to 3,000m above sea level	

Classification according to Paragraph 6 of the General Specifications for the Safety of Medical Electrical Equipment according to IEC 60601-1 / ANSI/AAMI ES 60601-1


ELECTRIC MICROMOTOR IM3 EM-12 L GA

The micromotor is used as described in the manual supplied by iM3 or visit www.im3vet.com

The Micromotor is intended to be used with contra-angles handpieces with couplings specified in ISO 3964, type 2 or type 3.

Lift up the micromotor and control it with the foot control.

To adjust the unit for intended use please see operation section of this manual which describes how to preset and adjust the motor.

To avoid necrosis, it is generally recommended to work with a 'wet' spray that gives at least 50ml spray per minute measured with the turbine handpiece fitted.

AIR INSTRUMENT

The hand piece connection is a 4 hole Midwest connection. iM3 range of recommended handpieces can be found on www.im3vet.com, manuals for each handpiece and recommended use can be found on www.im3vet.com.

To activate the high or low speed handpiece, lift up the handpiece and then press down on the foot pedal to operate.

The water coolant spray for the handpiece can be adjusted by using the blue coolant adjustment knob marked with the handpiece icon on the side of the dental unit as pictured below.

LOW & HIGH SPEED HANDPIECE PRESSURE

Handpiece operating pressures are preset in the factory, when pressing the foot pedal/operating the Low Speed motor the PSI should read 35psi-40psi and the Highspeed handpiece should read 40psi.

During your annual service your pressures will be checked by an iM3 Authorised Engineer, if the pressure needs altering outside this contact iM3 or your dealer www.im3vet.com

To remove the low speed attachment from the air-motor (low speed), pull the attachment upward away from the motor. To remove a high speed handpiece or a low speed air motor from the handpiece cable unscrew the plastic handpiece cable connection anticlockwise.

During use air will be forced back down the handpiece cabling, the air, is filtered before going into the atmosphere. The filter is located in the connection between the handpiece cable and the dental units and can be accessed by unscrewing the outer plastic as pictured below. This should be replaced at service or when oil is visible on the tubing below the connection.

To replace the cable twist to align the holes and then push upward.


Fig.15 - To remove the low speed attachment from the air-motor


Fig.16 - To remove/replace the handpiece cable from dental unit


Fig.17 - To replace the cable twist to align the holes and then push upward

ULTRASONIC SCALER

The Ultra Sonic Scaler is used as described by supplier, in the sub manual

Option 1 with the display attached.

Lift the scaler up and use the foot pedal to activate.

Use the plus and Minus keys on the digital display to increase and decrease the power of the scaler as per the recommended settings found on the Tip guide supplied with the Tip.

To adjust the water use the blue coolant adjustment knob with scaler icon (see figure 19)

 Please refer to the tip manual which accompanies the unit for the power settings for each tip.

Option 2

If the dental machine does not have a digital display use the scaler power adjustment, knob with the scaler icon to increase or decrease the ultrasonic scaler power (see figure 19)

To adjust the appropriate power for the scaling tip and the procedure please see **Operation section** of this manual which describes how to preset and adjust the scaler.

Two different Ultrasonic Scalers that can be attached to the Evolution dental unit. One option uses P6 style scaling Tips and the other option Ultra style scaling tips. The choice should be made at purchase. The ultrasonic scaler handpiece uses the symbols pictured (see figure. 20) , to indicate which tip should be used with the handpieces, tips selection can be found at www.im3vet.com.

Fig.18 - Ultrasonic scaler


Ultra style connection

P6 style connection

Fig. 20 - Ultrasonic Scaler different Tip types for handpieces

Fig.19 - Side panel controls

OVERVIEW OF DISPLAY

Factory Settings

- P** Program 1 Display shows P1: 30
- Program 2 Display shows P2: 20
- Program 3 Display shows P3: 15

- +** PLUS key: Increase power, maximum 40, possible during use
- MINUS key: Reduce power, minimum 1, possible during use

By holding PLUS/MINUS key the values are continuously increased / decreased.


Fig. 21 - Display Screen Ultrasonic Scaler

CHANGING / SAVING THE POWER ON THE ULTRASONIC SCALER

i The power can be set at P1, P2, P3 for a minimum of 1 to maximum 40

P Step 1: Press the program key (P1,P2 or P3).

+ Step 2: Increase power

- Step 3: Reduce power

P Step 4: Press program key P (P1, P2 or P3) for approximately 2 seconds to save.

i The Set values and an acoustic signal is emitted to confirm that all the settings have been saved.

ERROR MESSAGES FOR SCALER IF USING THE DISPLAY

DISPLAY	DESCRIPTION OF ERROR	SOLUTIONS
Error 1	Overheating / overloading of electronics	> ❶ Disconnect device from power supply > ❷ Wait 5 minutes and allow the system to cool > ❸ Restart the function
Error 2	Pedal of the foot control pressed during switching on	> Do not press the pedal of the foot control
Error 4	Display keys pressed during switching on	> Do not press the display keys
Error 5	Running time limiter as a result of 15 minutes of continuous operation	> Check the pedal of the foot control > Do not press the pedal of the foot control any more [release completely]
Error 6	Error with »Piezo Scaler«	> ❶ Check that the Piezo Scaler is correctly attached to the coupling > ❷ Replace the tips
Reboot		> Switch the power supply off and on again > Restart system
Error	e.g.: 05 6303	> Immediately contact an authorised iM3 service partner(see pg6)

SUCTION

To activate the suction, lift the handpiece out of the holder.
 To turn off, place the handpiece back into the holder.

The amount of suction can be adjusted on the suction handpiece (see figure 23).

If bleach or any bleach based product is used during treatment, then the suction tube should removed after treatment and cleaned fully. The tube should be then left to dry, hang straight not in a loupe or on the unit.

Fig.22 - Suction Instrument


Fig.23 - Suction can be adjusted on the suction handpiece


Fig.24 - Handpiece and Suction tubing can be disconnected for the purpose of cleaning

Fig.25 -The Collection bottle between Water and CLS Bottles

When Emptying the collection bottle, do not pour liquid over the Purple filter as this will stop the suction working and the filter will need to be replaced if this happen, Order code: P7957

DISTILLED WATER BOTTLE & CLS BOTTLE


Fig. 26 - Water & CLS Bottles


Fig. 27- Side of Unit controls

-  Your Evolution unit is supplied with a sterilising straw and a bottle of CLS, please follow the instructions for use before operating the dental unit. Watch the unboxing video supplied for more information on CLS and the straw, or visit www.im3vet.com for further information.
-   The water bottle is marked with a water droplet and the CLS bottle is marked with the words CLS. The CLS bottle is UV protected and a normal water bottle should not be used for CLS as it is light sensitive.
-  The water bottles will automatically pressurise when the bottles are attached to the dental unit and depressurise when they are detached from the unit.
-  Only tighten bottles hand tight, over tightening the bottles will reduce the life of the bottles and damage the bottle holder thread.
-  Only distilled/ionised water should be mixed with the CLS, if normal/tap water is used in conjunction with CLS then the unit can be blocked. Always follow the instructions on the CLS bottle for dilution amount.
-  Only Distilled/ionised water should be used in the water bottle, if normal/tap water is used in conjunction with the straw then the unit can be blocked.
-  If CLS is not used in the CLS bottle then the filter should be detached and a iM3 straw should be attached. If no straw is used then the machine will not be kept clean.

BUR REPLACEMENT - iPUSH TOOL

Bur Replacement - iPush Tool (iM3 ordering code L7050)

1. Place the yellow iPush onto the LED Advantage handpiece.
2. Depress the lever with your thumb and insert the bur all the way until it stops.
3. Remove the iPush before using the Handpiece.
4. To change the Bur repeat step 2. and remove the Bur. To increase the life of the Handpiece, Oil between each use and replace the Bur between each patient, (see Handpiece manual).


Fig.28 - Bur Replacement diagram

Handpieces included with your dental unit will depend on the model purchased. Visit www.im3vet.com to see which handpieces are included with your unit.

IM3 ADVANTAGE LOW SPEED HANDPIECE

The low speed handpiece features include a light weight balanced profile and smooth swivel action. Its vane type power system is quiet and vibration free. The motor accommodates both Doriot or ISO heads and its powerful torque is suitable for all low speed procedures. The Low speed Motor operates at 20,000 rpm with delivery pressure of 40 psi.

Operation Recommendations

- AIR PRESSURE** The dental unit air pressure should be adjusted to be 40 psi
- LUBRICATION** To assure optimal performance and life, the handpiece should be lubricated daily with MD-30 oil lubricant.
- AIR SUPPLY** Moist or dirty air from the compressor can cause damage. Check the regulator filter regularly. If dirt or moisture is present, replace the regulator filter and bleed the system, to remove water

Installation Instructions

1. Connect the low speed handpiece to the grey handpiece line
2. Fit the 4:1 Advantage straight nose cone. Attach the Blue R&R oscillating prophy head by simply pushing it on
3. Depress the foot control fully to operate the low speed handpiece

Fig.29 - Low Speed Handpiece

Polishing should be carried out at 5,000 rpm or less. This can be accomplished by running the motor at full revs (20,000 Rpm) and use a 4:1 reduction head (L6250), to achieve 5,000 revs, maintaining torque.

iM3 Low Speed Handpieces Attachment

(1) 4:1 Advantage Straight Nose Cone (SNC) Attachment L6250

The reducing head is designed for polishing (prophy). The speed of the polishing cup should be maintained at less than 5,000 rpm. The reducing attachment allows the motor speed and torque to be maintained while reducing the polishing speed. iM3 disposable oscillating prophy angles are fitted to the SNC. They won't trap muzzle hair, or generate excessive heat. (iM3 R&R prophy cup Code L7466)



Fig.30 - Advantage Straight Nose Cone attachment

Fig.31 - Advantage Straight Nose Cone

Straight Nose Cone Attachment Blue Band For HP burs and used with the iM3 rodent kit. This includes a soft tissue protector with diamond bur and diamond disc.

iM3 HIGH SPEED HANDPIECE

The Advantage LED High Speed (400,000 rpm) Handpiece and Swivel (L7300) with inbuilt generator creates its own bright LED light. With a colour temperature of 5,500 K the Advantage produces a radiation spectrum that corresponds to daylight. The high-colour rendering index (CRI 90) of the LED enables relaxed work thanks to daylight-quality light and provides exceptional reproduction of natural colours and clear colour contrasts in the mouth. The Advantage minimizes fatigue thanks to the easy 360 degree rotation of the handpiece on the swivel coupling, the low weight and the slender design. The unique handpiece design also allows for component replacement on site by the user under guidance from iM3

All handpieces should be lubricated prior to use with iM3 MD-30

Cat. Code MD30

Fig.32 - iM3 MD-30 Oil can be used to lubricate your stone.


Fig.33 - High Speed Handpiece

iM3 BOTTLE SYSTEM

INSTALLATION & OPERATING INSTRUCTIONS

Fig.34 - Units CLS and Water bottles

CONTENTS

- › Installation & Operating Instructions 31
- › 1. Safety 32
- › 2. Accessories 34
- › 3. Technical data 35
- › 4. Function 36
- › 5. Handling 37
- › 6. Disinfection and cleaning 37
- › 7. Reprocessing 39
- › 8. Maintenance 42
- › 9. Tips for operators and service technicians 43
- › 10. Manufacturer of Bottles 43

1. SAFETY

iM3 has developed and designed the system in such a way that dangers are effectively ruled out if the unit is used in accordance with the intended use.

Despite this, the following residual risks can remain:

- » Personal injury due to incorrect use/misuse
- » Personal injury due to mechanical effects
- » Personal injury due to lack of hygiene, e.g. infection

1.1 Intended Purpose

The iM3 Bottle System is used to supply water to dental treatment units and in particular to the dental instruments connected to them.

1.2 Intended Use

Fill the iM3 bottle with a suitable solution. This can be a disinfection solution to prevent algae formation and to disinfect the waterways with drinking water or a suitable sterilising solution.

The system must be connected to a compressed air line as well as to the line supplying the treatment unit with water (see "6 Mounting"). The system may only be used with the supplied safety pressure regulator. The integrated stop valve enables the bottle to be screwed in and out while the system is switched on.

1.3 Improper Use

Any other usage or usage beyond this scope is deemed to be improper. The manufacturer accepts no liability for damage resulting from improper usage. In such cases, the user/operator will bear the sole risk. Do not use any liquids or solids that are unsuitable for the iM3 Bottle System.

1.4 Systems, connection with other devices

Additional devices connected with medical devices must be proven to conform with their corresponding IEC or ISO standards. All configurations must continue to comply with the standard requirements for medical systems (see IEC 60601-1).

Whoever connects additional devices to medical devices automatically becomes the system configurator and is responsible for ensuring that the system corresponds with the standard requirements for systems. Local laws take priority over the requirements outlined above.

1.5 General safety notes

- » Always comply with the specifications of all guidelines, laws, and other rules and regulations applicable at the site of operation for the operation of this unit.
- » Check the function, condition and fill level of the unit prior to every use.
- » Do not convert or modify the unit.
- » Comply with the specifications of the Installation and Operating Instructions.
- » The Installation and Operating Instructions must be accessible to all operators of the unit at all times.
- » When filling the bottle, watch out for spilled liquids - danger of slipping.

1.6 Combining devices safely

Where applicable, the requirements for medical products have been taken into account in the development and construction of the device. As a result, this device is suitable for installation within medical supply equipment.

- » Where this device is integrated in other medical supply equipment, the requirements of European Union Medical Device Regulation (EU) 2017/745 and the relevant standards must be observed.

1.7 Specialist personnel

Handling

Unit operating personnel must ensure safe and correct handling based on their training and knowledge.

- » Instruct or have every operator instructed in handling the unit.

The following groups are not permitted to operate or use a commercially operated unit:

- » People without the necessary experience and knowledge
- » People with reduced physical, sensory or mental capabilities
- » Children

Installation and repairs

- » All installation, resetting, alteration, extension and repair work must be carried out either by iM3 personnel or by a suitably qualified person approved by iM3.

1.8 Essential performance characteristics

The unit does not have any essential performance characteristics as set out in IEC 60601-1 section 4.3. The unit complies with the requirements according to DIN EN ISO 60601-1: 2022 Cap. 9.7 and DIN EN ISO 7494-2:2023-03.

1.9 Duty to report serious incidents

i The operator/patient is required to report any serious incident that occurs in connection with the device to the manufacturer and to the competent authority of the Member State in which the operator and/or patient is established/resident.

1.10 Only use original parts

- » Only use accessories designated or approved by iM3.
- » Only use only original wear parts and replacement parts.

1.11 Storage and transport

The original packaging provides optimum protection for the unit during transportation.

- i** iM3 will not accept any responsibility or liability for damage occurring during transportation due to the use of incorrect packaging, even where the unit is still under guarantee.
- » Only transport the unit in its original packaging.
 - » Keep the packing materials out of the reach of children.

i The bottles are to be stored UV-protected when not in use.

1.12 Disposal

- i** The unit may be contaminated. Inform the waste management company so that they can take all necessary safety steps.
- » Prepare accessory parts before disposal, then dispose of according to local and national regulations.
 - » If you have any questions about the correct disposal of parts, please contact your dental trade supplier.

Fig.35 - iM3 Evolve Water bottle

2. ACCESSORIES

2.1 Bottle Parts

NOTICE

Device malfunction or damage due to incorrect/missing accessories
 Guarantee claims may become invalid as a result. The bottles can burst.

- » Only use iM3 bottles.
- » Only operate with a safety pressure regulator 2.0 bar.
- » Exchange bottle after one year or 50 automatic treatments.
- » Damaged, scratched or opaque bottles must be replaced immediately.

i An overview of iM3 bottles is available on the iM3 online shop: www.im3vet.com

Water Bottle with iM3 logo, lasered **Order code: BT001**

CLS Bottle with iM3 logo, printed **Order code: BT002**

i Other inscriptions on the bottle are available on request.

i The safety pressure regulator for fixed hoses requires a hose with an outer diameter of 4 mm.

2.2 Wear and replacement parts

The following parts are subject to wear and tear, and should be replaced at regular intervals (see also **10 Maintenance**):

Riser pipe with fine filter: **Order code: RM012**

Straw filter **Order code: PS001**

i The safety pressure regulator must be replaced at intervals of 4 years.

Contact iM3 support for spare parts:

i EU & UK: support@im3dental.com
 US: info@im3usa.com
 AU: sales@im3vet.com.au

3. TECHNICAL DATA

Medical Device Class		
iM3 Bottle System		Class 1
General technical data		
Dimensions (W x H x D)	cm	approx. 9 x 30 x 10
Weight	kg	approx. 0.5
Hose diameter		
PUR air hose blue (ID x OD x L)	mm	1.5 x 3.45 x 500
PUR water hose green (ID x OD x L)	mm	1.7 x 3.2 x 500
Water flow rate		
Standard variants	ml/min	approx. 760
Cup filler	ml/min	approx. 2155
Cup filler hose dimensions		
PUR air hose blue (ID x OD x L)	mm	1.5 x 3.45 x 500
PUR water hose green (ID x OD x L)	mm	3.0 x 4.5 x 500
Ambient conditions during operation		
Temperature	°C	+3 to +50
Relative humidity	%	5 to 95
Minimum connection pressure – air	bar	2
Maximum connection pressure – air	bar	6
Ambient conditions during storage and transport		
Temperature	°C	-20 to +70
Relative humidity	%	5 to 95
Atmospheric pressure	hPa	800 to 1200

3.1 Type plate

The type plate with the UDI code of the iM3 Bottle System can be found on the back of the bottle receptacle.

3.2 Evaluation of conformity

This unit has been subjected to a conformity assessment procedure in accordance with the current relevant European Union guidelines. This unit conforms to all relevant requirements.

4. FUNCTION

The unit is used to supply water to dental treatment units and in particular to the dental instruments connected to them.

The system must only be connected to a iM3 bottle.

Another function is to disinfect the waterways in the dental instrument, including the iM3 Bottle System, and/or to remove and prevent algae growth, using a disinfectant solution in the bottle. The bottle is filled with the appropriate medium (drinking water or disinfection solution). The system is affixed to a compressed air line and to the line that provides the treatment unit with water and then fastened appropriately (see 6 Mounting). Then the system is operational.

When unscrewing the bottle, the attached compressed air line is blocked by a stop valve, which is re-opened when screwing in the bottle. This makes it easy to change bottles or to refill the bottle while the unit is operational.

Fig.36 - Installing Units bottles

A	iM3 bottle, screwed in
B	iM3 bottle, unscrewed

5. HANDLING

5.1 Initial use

- » Clean the iM3 Bottle System before initial use (see "8 Disinfection and cleaning").

5.2 Operating tips

- » Insert the riser pipe into the bottle.
- » Place the bottle in the bottle receptacle.
- » Turn anti-clockwise approx. 1/2 turn and tighten gently. The pressure is mounted automatically.

Fig.37 - Mounting the Units bottles

- » Insert riser pipe in the bottle [A]
- » Screw in the bottle [B]

6. DISINFECTION AND CLEANING

WARNING

Health hazard due to contaminated unit
Contaminated units can cause infections.

» Clean and disinfect the iM3 Bottle System before working on the unit.

» Wear suitable personal protective equipment when working (e.g. impermeable gloves, protective goggles and mouth and nose protection).

NOTICE

Device malfunction or damage due to incorrect cleaning and disinfecting agents

Guarantee claims may become invalid as a result.

- » Do not use any foaming agents, e.g. household cleaning agents or instrument disinfectants.
- » Do not use abrasive cleaners.
- » Do not use agents containing chlorine.
- » Do not use any solvents like acetone.

NOTICE

Device malfunction or damage due to damaged/worn bottles

Guarantee claims may become invalid as a result.

The bottles can burst.

Bottles are subject to operational wear and tear and have a limited service life*.

- » Exchange bottle after one year or 50 automatic treatments.
- » Damaged, scratched or opaque bottles must be replaced immediately.

* The expiry date is lasered on the bottom of the bottle, along with the corresponding pictogram.

Fig.38 -Bottles Expiry date on bottle

iM3 recommends:

- ✓ For continual disinfection of procedural water and water supply lines: **Alpron* solution**
- ✓ For cleaning and disinfection of waterways: **Bilpron* solution**
- ✓ For cleaning and disinfection of bottles: **BC-San 100* solution**
- ✓ For cleaning and disinfection of external surfaces, bottles and threads: **FD 366 sensitive top wipes Disinfection of sensitive surfaces and FD 366 sensitive Disinfection of sensitive surfaces, ready-to-use solution****

Only these products have been tested by iM3.

* Alpron, Bilpron and BC-San 100 are products from Alpron Medical GmbH.

** FD 366 sensitive top wipes and FD 366 sensitive Disinfection of sensitive surfaces, ready-to-use solution, are products from DÜRR DENTAL SE.

6.1 Procedural water disinfection

Drinking water can be mixed with Alpron for use on patients. This ensures consistent water quality.

- » Detailed information can be found in the Alpron operating instructions.

6.2 Daily after the end of treatment

- i** The bottle contents must be changed every day, otherwise there is a risk of contamination.

The following are required for cleaning:

- » Non-foaming disinfectant/cleaning agent that is compatible with the materials.
- » Clean the bottle with hot water (≤ 90 °C) and a cleaning brush.

Cleaning and disinfecting surfaces

- » Regularly clean and disinfect surface and thread of the bottle and bottle receptacle, e.g. with FD 366 sensitive top wipes.

6.3 Once or twice a week/break-time preservation

The following are required for cleaning:

- » Non-foaming disinfectant/cleaning agent that is compatible with the materials.
- » Rinse the bottle with hot water.
- » Clean the bottle with a cleaning brush or in the washer-disinfector.
- » Fill the bottle with Bilpron according to the manufacturer's instructions and flush all the lines. Detailed information can be found in the Bilpron operating instructions.
- » Fill the empty bottle with BC-San 100 according to the manufacturer's instructions and close with a screw cap, if available.

- i** A bottle that is filled with BC-San 100 solution must not be connected.

Detailed information can be found in the BC-San 100 operating instructions.

7. REPROCESSING

7.1 Risk analysis and categorisation

A risk analysis and categorisation of medical products often used in dentistry must be performed before their reprocessing by the operator. Comply with all national directives, standards and specifications such as, for example, the "Recommendations from the Commission for Hospital Hygiene and Infection Prevention".

Accessories of the medical device are also subject to reprocessing.

Classification recommendation given proper use of the product: semi-critical B.

Semi-critical medical device

A medical product that comes into contact with mucous membranes or pathologically affected skin. The operator is responsible for correct classification of the medical products, defining the reprocessing steps and performing the reprocessing.

7.2 Preparation process in accordance with EN ISO 17664

- » The reprocessing procedure after each treatment is carried out according to the reprocessing procedure established by EN ISO 17664.

Important information!

The reprocessing notes in accordance with EN ISO 17664 have been independently tested by iM3 for the preparation of the device and its components for their reuse.

The person conducting the reprocessing is responsible for ensuring the reprocessing performed using the equipment, materials and personnel achieves the desired results. This requires validation and routine monitoring of the reprocessing process. Any deviation from the procedures described following by preparation staff is solely the responsibility of the personnel in terms of effectiveness and the possible adverse consequences.

Frequent reprocessing has little effect on the device components. The end of the product life cycle is especially influenced by the amount of wear and tear or damage resulting from its use. The use of soiled, contaminated and damaged components is at the sole responsibility of the person performing the reprocessing and the operator.

The reprocessing method was validated as follows:

- » Pre-cleaning
 - » FD 366 sensitive top wipes
 - » Cleaning brush
 - » Manual cleaning
 - » FD 366 sensitive top wipes
 - » FD 366 sensitive ready-to-use solution
 - » Cleaning brush
- » Manual disinfection
 - » FD 366 sensitive top wipes
 - » Bilpron
 - » BC-San 100 (only for the bottle)
- » Automatic cleaning and disinfection
 - » was performed in accordance with EN ISO 15883 with tested efficacy.
 - » Cleaning agent: Neodisher MediClean Forte
 - » Washer-disinfector: PG 8535 (Miele)
 - » Programmes: "Cleaning without neutralisation" and "THERMAL DESotro"
 - » Cleaning/disinfection temperature: 90 °C
 - » Cleaning/disinfection time: 5 minutes
- » Cleaning brush
 - » Bottle brush with silicone bristles
 - » Bristle material: silicone
 - » Brush head: Ø 60 mm
 - » Temperature resistance up to 300 °C
 - » Brush length: 350 mm

Example:

Culinaris bottle brush, sterilisable, anthracite REF 227721

General information

- » Comply with all national directives, standards and specifications for the cleaning, disinfection and sterilisation of medical products as well as the specific specifications for dental practices and clinics.
- » When selecting the cleaning and disinfecting agents to be used, observe the specified information (see “Manual cleaning, intermediate rinsing, disinfection, final rinse, drying” and “Automatic cleaning, intermediate rinsing, disinfection, final rinse, drying”).
- » Comply with the concentrations, temperatures, residence times and post-rinsing specifications issued by the manufacturer of the cleaning and disinfectant agent.
- » Only use cleaning agents that are non-protein-fixing, aldehyde-free and display material compatibility with the product.
- » Only use disinfectants that are aldehyde-free and display material compatibility with the product.
- » Do not use any rinse aid (danger of toxic residue on the components).
- » Only use freshly-produced solutions.
- » Only use distilled or deionised water (not sterile water) with a low bacterial count (at least drinking water quality) that is free from facultative pathogenic microorganisms (e.g. legionella bacteria).
- » Use clean, dry, oil and particle-free compressed air.
- » Do not exceed temperatures of 93 °C.
- » All devices used (e.g. ultrasonic bath, cleaning and disinfection device (washer-disinfector) must be regularly serviced and inspected.

7.3 Preparation at the operating location

Wear protective gloves

Wear protective goggles.

Use mouth and nose protection

Use protective clothing

Fig.39 - Riser piper

WARNING

**Risk of infection from contaminated products.
Risk of cross contamination**

- » **Prepare the product correctly and promptly before its first use.**

- » Unscrew the bottle from the bottle receptacle [A].
- » Remove the riser pipe from the bottle receptacle [B].
- » Thoroughly wipe down the exterior surfaces of all components with cleaning wipes to remove coarse organic and inorganic soiling:
 - » one cleaning wipe for the small components, e.g. the individual parts of the riser pipe and
 - » two cleaning wipes for larger components, e.g. the bottle.
- » Note the action time of the cleaning agent.
- » Protect the unit from contamination when transporting it from the treatment chair to the reprocessing location.

7.4 Manual cleaning, intermediate rinsing, disinfection, final rinse, drying

A disinfectant or combined cleaning and disinfectant agent is required for manual disinfection. It must have the following properties:

- » Proven effectiveness (VAH or European standards) For more information, see "General information".

Cleaning

- » Thoroughly wipe riser pipe and bottle receptacle as well as their threads with cleaning and disinfection wipes (non-protein-fixing/aldehyde-free, see "General information"), until all contamination has been removed.
- » Wipe bottle with cleaning and disinfection wipes, clean inner surfaces and thread with cleaning brush and cleaning and disinfecting solution, to remove all dirt.
- » Comply with the reaction times of the cleaning and disinfecting agents (see "General information").
- » If you notice any further contamination, brush all surfaces completely with a cleaning brush and cleaning and disinfecting agent.

Intermediate rinsing

After the action time prescribed by the manufacturer:

- » Rinse all components under running water (temperature < 35 °C) for at least 1 minute.
- » Disinfect
- » Wipe and disinfect the riser pipe and bottle receptacle, as well as their threads, with cleaning and disinfecting agent (non-protein-fixing/aldehyde-free, see "General information").
- » Wipe and disinfect the exterior surface of the bottle with cleaning and disinfecting wipes.
- » Fill the bottle with Bilpron, connect it to the iM3 Bottle System and flush all the lines. Detailed information can be found in the Bilpron operating instructions.
- » Fill empty bottles that are not connected to the unit with BC-San 100 then close. Detailed information is available in the BC-San 100 operating instructions.
- » Note the action time for the disinfectant.

Final rinse

After the action time prescribed by the manufacturer:

- » Rinse all components under running water (temperature < 35 °C) for at least 1 minute.

Drying

- » If necessary, re-dry components in a clean location using a hygienic, lint-free cloth.

- » Dry the components with compressed air in a clean location.
- » Store empty bottles with the opening facing downwards – observe drying time of 4 hours.

7.5 Automatic cleaning, intermediate rinsing, disinfection, final rinse, drying

Selection of the washer-disinfector

Automatic cleaning and disinfection requires a washer-disinfector with the following properties and validated processes:

- » Corresponds to and tested in accordance with EN ISO 15883
- » Inspected programme for thermal disinfection (A0 value ≥ 3000 or min. 5 minutes at 93 °C)
- » Programme is suitable for the components and provides sufficient rinsing cycles (see also "General information").

Selection of the machine cleaning agents and disinfectants

The following properties are required:

- » Material compatibility with the product
- » Corresponds with the manufacturer's specifications of the CD

For further information, see: "General information").

Automatic cleaning and disinfecting

When arranging the parts in the washer-disinfector, make sure there are no areas missed by rinsing.

Observe the drying time for the bottle of min. 4 hours.

- » Attach bottle to suitable holders in the washer-disinfector.

7.6 Check for function

- » After the end of the cleaning and disinfection cycle, check the components for any residual soiling and residual moisture. If necessary, repeat the cycle after the drying time.
- » If necessary, replace any damaged parts, see "3.3 Accessories".
- » Insert the riser pipe in the bottle receptacle, screw in the bottle and check the system for leaks.
- » The components should be packaged as soon as possible after drying and checking.

8. MAINTENANCE

All maintenance work must be performed by a qualified expert or by one of our service technicians.

WARNING

Health hazard due to contaminated unit
Contaminated units can cause infections.

» Clean and disinfect the iM3 Bottle System before working on the unit.

» Wear suitable personal protective equipment when working (e.g. impermeable gloves, protective goggles and mouth and nose protection).

8.1 Replacing the riser pipe


Fig. 40 - Replacing the riser pipe

- » Remove the riser pipe [A].
- » Insert a new riser pipe [B].

Riser pipe (1 piece) **Order code: RM012**

9. TIPS FOR OPERATORS AND SERVICE TECHNICIANS

 **Any repairs exceeding usual routine maintenance must be carried out by suitably qualified personnel.**

 **WARNING**
Infection due to contaminated unit
 Contaminated units can cause infections.

- » Clean and disinfect the iM3 Bottle System before working on the unit.
- » Wear suitable personal protective equipment when working (e.g. impermeable gloves, protective goggles and mouth and nose protection).

ERROR	POSSIBLE CAUSE	SOLUTION
Low water flow rate	Fine filter blocked in the riser pipe	Replace the fine filter.
	Blockage in the riser pipe	Remove and clean riser pipe.
	Blockage in the PUR hose	Remove and clean the PUR hose, replace if necessary.
	Kink in the PUR hose	Reposition the PUR hose.
No water flow	Blockage in the riser pipe	Remove and replace riser pipe.
	Blockage in the PUR hose	Remove and replace the PUR hose.
	No pressure	Switch on the compressed air supply. Screw in the bottle as far as it will go, to open the stop valve.
Fluids leakage from the bottle system	PUR hose is not attached correctly	Reconnect the PUR hose with a hose safety device.
Air blows off audibly	Bottle not completely screwed in	Turn bottle to screw in completely.
	Safety valve blows off	Pull the ring of the safety valve to vent.

10. MANUFACTURER OF BOTTLES

 **Manufactured and serviced by:**

iM3 Dental Limited
 Unit 29, Duleek Business Park,
 Duleek, Co. Meath,
 A92 N72W, Ireland.
 Phone: +353 (0) 1 691 1277
 Email: sales@im3dental.com |
 support@im3dental.com
 Website: www.im3vet.com

iM3 Inc.
 12414 NE 95th St., Vancouver, WA,
 98682, USA
 Phone: +1 800 664 6348
 Fax: +1 360 254 2940
 Technical Support: +1 800 664 6348
 Email: info@im3usa.com
 Website: www.im3vet.com

iM3 Pty Ltd.,
 21 Chaplin Drive, Lane Cove, NSW
 2066, Australia.
 Phone: +61 2 9420 5766
 Fax: +61 2 9420 5677
 Email: sales@im3vet.com
 Website: www.im3vet.com.au

4. INFECTION CONTROL

CROSS-CONTAMINATION

The following procedures are recommended to minimize the risk of cross-contamination

If the described cleaning, disinfection and sterilization methods are not followed careful: The safety of operators and patients may be compromised, the service life of the unit will be compromised and iM3 has no responsibility for the products correct functioning and safety.

Before using your Evolution Dental Unit for the first time it must be cleaned, disinfected and sterilized as described below:

- » Before using a Thero disinfector and autoclave, read the instructions for use of the devices and be aware of the warnings provided by the manufacturers.
- » For any instruments other than Xo products, always follow the cleaning, disinfection and sterilization instructions included by the manufacturer.
- » Pack items to be sterilized in sterilization packages that meet the following requirements:
 - The Sterilization package must meet the applicable standards in respect of quality and use, it must be suitable for sterilization method
 - The Sterilization package must not be under tension
- » Store sterile goods dust-free and dry.

The following definitions apply:

Cleaning, Disinfection & Sterilization definitions

PROCEDURE	PURPOSE	METHOD
Cleaning	Removal of visible blemishes, spots and stains etc.	<ul style="list-style-type: none"> • Physical/chemical use of soap or detergents, will not significantly reduce the number of pathogens
Disinfection	Significantly reduce the number of pathogenic micro-organisms	<ul style="list-style-type: none"> • Application of chemical disinfection to surfaces • Application of thermos-disinfection according to ISO 15883-1:90°C /194°F, 5 minutes.
Sterilization	Elimination/destruction of all living pathogenic micro-organisms	<ul style="list-style-type: none"> • Steam sterilization (steam sterilizer according to EN 13060, Type B) • 134°C/273°F for 3 minutes or • 121°C/250°F for 15 minutes

Please note that some instruments and accessories are not designed for thermo-disinfectors or autoclaves!
 Please note that autoclaving and thermo-disinfection process wear down the materials and may cause colour change and/or shorten the life time.

CLEANING, DISINFECTION & STERILIZATION PROCEDURES

Except for mild soap or detergents (e.g. liquid soap) please only use detergents and disinfectants approved by iM3, a full list can be found on www.iM3vet.com

Always follow the instructions provided in the instructions of the product for concentration and contact time.

Please note: Do not use any other disinfectants than the ones listed on the iM3 website. Doing so may damage the unit. Failure to comply with these precautions may affect the Evolution product warranty. If in doubt contact iM3.

GENERAL DISINFECTION OF UNIT SURFACES, HANDLES AND SKIRT

Disinfect the surfaces of the unit using a lint-free cloth with recommended disinfection. Look out for disinfection liquid dripping from the cloth.

When disinfecting surfaces – do not use an excessive amount of disinfectant and make sure that the disinfectant has evaporated before putting silicon mat and handpieces back into place.

INFECTION CONTROL OF THE SILICON MAT

The silicone mat can be cleaned and disinfected in the autoclave or thermo-disinfector


Fig.43 - Evolve Veterinary Dental Unit Disinfect all elements

DISINFECTION OF INSTRUMENT HOSES

Do not wash instrument hoses in the autoclave or theros-disinfector

Display:

The device is not approved for mechanical cleaning (thermo washer disinfector) and sterilization (Autoclave). Do not immerse in water or clean it under running water

Pre-disinfection:

If heavily soiled, clean first with disinfectant cloths.

Only use disinfectants that have no protein fixing effects

Use only disinfectants which do not contain chlorine and

INFECTION CONTROL OF OTHER INSTRUMENTS

For turbines, ultra sonic scalers and other instruments please refer to their user manuals supplied with the units.

DISINFECTION OF UNIT WATER LINES

The iM3 Straw should be used at all times on the water line of your dental units. As standard your unit will have come with a Straw, please follow the instructions carefully on how to use this item.


It is a violation of federal/local law to use this product in a manner inconsistent with its labeling.

INITIAL INSTALLATION:

Unscrew the water bottle on your dental unit.

- » Fill bottle with distilled quality water.
- » Fill the bottle, leaving an air space of at least 1" from the top of the bottle (do not over fill the bottle - adequate airflow is necessary to ensure proper pressurization).


The iM3 Straw includes a Luer Lock.

Proper installation of the luer lock is shown in the diagrams. It is important to cut the existing pickup tube to the proper length so that the bottom of the iM3 Straw is at least a 1/4" from the bottom of the bottle.

INITIAL TREATMENT: iM3 Straws are engineered to automatically produce a "Shock" treatment after initial connection and treatment of dental water. Slowly purge the waterlines attached to the dental unit into a cup until you see the light blue colour. The colorant is an indicator the Shock product is in the line(s). Allow lines to remain unused for a minimum of 12 hours. For best results, allow Shock to remain in the line(s) over the weekend. Prior to use of dental unit, flush dental lines until blue colour disappears (approximately 4 ounces).

MAINTENANCE TREATMENT: Fill water bottle with distilled water. The dental water is automatically treated as it runs through the Straw Cartridge. Replace the Straw Cartridge per the Straw Capacity & Replacement Schedule, at least yearly.

HAZARDS TO HUMANS AND DOMESTIC ANIMALS - CAUTION:

In the event of a rupture of the cartridge, avoid contact with the silver ion exchange resin which may cause eye and skin irritation. Avoid direct contact with eyes, skin and clothing. Wash thoroughly with soap and water after handling and before eating, drinking, chewing gum, using tobacco or using the toilet. Remove and wash contaminated clothing before reuse. Place debris into a sealed plastic bag and discard in trash.

ENVIRONMENTAL HAZARDS: This pesticide is toxic to fish and aquatic invertebrates.

FIRST AID: If in eyes: Hold eye open and rinse slowly and gently with water for 15–20 minutes. Remove contact lenses (if present) after the first 5 minutes then continue rinsing. If on skin: Take off contaminated clothing. Rinse skin immediately with plenty of water for 15–20 minutes. Call a poison control centre or doctor for treatment advice.

Handpiece

Refer for more details to the handpiece operating manual.

which are certified by officially recognized institutes.
Note the manufacturers specification for the use of the disinfectants

Before resuming operation wait until completely dry

HIGH SPEED HANDPIECE

Handpieces should be thoroughly cleaned to remove any foreign particles, eg ; blood, saliva . Make sure that the turbine cartridge and head cavity are clean before autoclaving. If autoclaved with bloody particles in the turbine housing, autoclaving will cause particles to be scorched and damage to the turbine may result.

The exterior of the handpiece may be cleaned with any good surgical disinfectant. Be sure to rinse and dry thoroughly. Lubricate handpiece turbine with MD-30 Lubricant. Do not submerge the handpiece in liquid. Refer to handpiece operating manual.

LOW SPEED HANDPIECE

Prepare handpiece and attachments for sterilization by carefully cleaning the exterior surface. Thoroughly scrub, rinse and clean away any residual solution and particles. Remove moisture with a towel or use air from the 3 way air water syringe.

Prior to sterilization, lubricate the handpiece using MD-30 Lubricant. Do not submerge the handpiece in liquid

HANDPIECES: Apply oil through the handpiece air inlet tube. As pictured on the iM3 MD-30 bottle

ATTACHMENT: Remove attachment from handpiece and oil through both end openings.

After lubricating items, assemble and operate for a few seconds to expel excess lubricant

1. Seal items in sterilization tubing bags. Make sure to remove bur from handpiece prior to sealing.
2. Place handpiece or attachments in sterilization trays of dry heat sterilizer or steam autoclave. Do not exceed temperatures of 275°F or 135°C. Sterilize according to manufacturer's instructions.
3. After sterilization cycle is completed, allow handpieces and attachments to cool down then lubricate handpiece and attachments thoroughly. Make sure to use a different container after sterilization cycle, to prevent re-contamination.

Lubricate handpieces using MD-30 Lubricant after cooling down, do not allow handpiece to sit in sealing bag overnight or for a long period. Immediately operate handpiece as described in the following step.

4. After lubricating items, insert bur into handpieces and operate for a few seconds to remove excess lubricant.

Do not autoclave the high speed handpiece with a bur in the chuck.

ULTRASONIC SCALER

For cleaning information read in full the sub-manual for the Scaler

Fig.44 - Handpiece Bridge

Fig.45 - Prior to sterilization lubricate the handpiece using iM3's MD-30

5. UNIT MAINTENANCE

TECHNICAL DATA UNIT

EMC Testing

DISPLAY

MODEL	D-2.0 E
Power supply:	5 V DC \pm 10%, from W&H controller ---V
Power consumption:	1.5 VA
Dimensions in mm (WxDxH):	155.5 x 29 x 122
Weight:	280 g

✓ AMBIENT CONDITIONS

Temperature during storage and transport:
-40°C to +70°C (-40°F to +158°F)

Humidity during storage and transport:
8% to 80% (relative), non-condensing

Temperature during operation:
+10°C to +35°C (+50°F to +95°F)

Humidity during operation:
15% to 80% (relative), non-condensing

Permitted ambient pressure: 70 – 106 kPa

Pollution level: 2

Overvoltage category: II

Altitude: up to 3,000 m above sea level

ULTRASONIC SCALER HANDPIECE

	PB-5 L, PB-5 LS, PB-5 LQ
Max. power output to the handpiece with load (ultrasonic) (W)	10
Frequency (ultrasonic) (kHz)	22–35
Minimum coolant supply volume (ml/min)	0*/20
Maximum coolant supply volume (ml/min)	50
Water pressure (bar)	1–6
Max. oscillating amplitude (Tip 1U) (mm)	0.2

* for tips where no coolant is used

Classification according to § 6 of the General Specifications for the Safety of Medical Electrical Equipment according to IEC 60601-1/ANSI/AAMI ES 60601-1

Type B applied part
(not suitable for intracardiac application)

TEMPERATURE INFORMATION

Temperature of the medical device at the operator side:
maximum 71°C (159.8°F)

Temperature of the medical device at the patient side
(front area of the medical device): **maximum 50°C (122°F)**

Temperature of the medical device at the patient side
(optical fibre): **maximum 48°C (118.4°F)**

Temperature of the working part (tip): **maximum 41°C (105.8°F)**

✓ AMBIENT CONDITIONS

Temperature during storage and transport:
-40°C to +70°C (-40°F to +158°F)

Humidity during storage and transport:
8% to 80% (relative), non-condensing

Temperature during operation:
+10°C to +35°C (+50°F to +95°F)

Humidity during operation:
15% to 80% (relative), non-condensing

Pollution level: 2

Overvoltage category: II

Altitude: up to 3,000 m above sea level

MICROMOTOR

Motor	EM-11 L	EM-12 L
Approved supply hose	VE-11	VE-10 / VE-11
Transmission instrument according to standard	ISO 3964	
Direction of rotation	forward/reverse	
Speed range	2,000 – 40,000 rpm	100 – 40,000 rpm
Maximum torque at the motor	3 Ncm	
Adjustment cooling air	6 – 8 NI/min	
Air coolant pressure* The air coolant pressure has to be higher than the water coolant pressure	0.5 – 3.0 bar	
Water coolant volume at (0,5 bar)	> 60 ml/min	
Water coolant pressure*	0.5 - 3.0 bar	

* Adjust the actual pressure with an attachment in place.

Supply hose	VE-10	VE-11
Approved electric motor	EM-12 L	EM-11 L / EM-12 L
Drive air respective cooling air at 250 kPa (2,5 bar)	> 8 NI/min	
Spray air at 250 kPa (2,5 bar)	> 8 NI/min	
Spray water at 200 kPa (2,0 bar)	> 200 ml/min	
Maximum pressure	400 kPa (4.0 bar)	

COMPRESSOR

Rate Voltage	AC 220V/50HZ
Power	<840W @ 700KPa
Rated Pressure	700KPa
Rate Flow L/min	>57 l/min
Ambient degrees C	5 – 40° C
Insulation Class	B
Thermal switch	Auto reset 145+/- 5 C
Net weight kg	<11.5kg
Restart pressure kPa	0 kPa
Tank Size	11 litres
Tank pressure	110 psi/7.5BAR maximum operating pressure

Ensure that the parts are not contaminated on disposal.

Follow your local and national laws, directives, standards and guidelines for disposal.

> Waste electrical equipment

> Packaging

6. UNIT ANNUAL SERVICE

EVOLVE DENTAL UNIT

iM3 products have been manufactured with great care by highly qualified technicians using a wide variety of tests and controls designed to ensure faultless operation.

To maintain a high standard of quality control, every dental system sold by iM3 is thoroughly tested and operated complete with the hand pieces and coolant before it is packaged for delivery. The quality of workmanship and unique design features used to build every iM3 are so good that our range of dental machines now carry a warranty which is setting new standards within the veterinary dental industry. Up to three years for our range of air-driven dental machines.

At iM3, our reputation and passion is your guarantee of advanced technology, unique features and quality workmanship. We design them, we build them and we service them. We know your iM3 workstation inside and out.

Technical & Support

Call us anytime during business hours on our Tech Support Line and we'll take care of all your needs or email us with your query and we'll get back to you promptly

For servicing quotation contact us at:

Email: support@im3dental.com

Phone (Ireland): +353 (0) 1 691 1277

Phone (UK): +44 (0) 1 423 224297

or message us via our Website online form:

<https://im3vet.eu/servicing-eu/>

(USA): 1 800 664 6348 / info@im3usa.com

(AUS): +61 2 9420 5766 / sales@im3vet.com

FREE Technical Support is included with every iM3 Dental Workstation

The relationship with our customers is of the utmost importance to us, to ensure that our products are performing to a standard we're proud to offer. The iM3 superior Technical Support Team comes FREE with every iM3 dental workstation or iM3 product.

Our Technical Support Team will continue to look after your dental system and all questions regarding handpieces, operational issues or maintenance long after you've taken delivery. We will help you to achieve the most out of your iM3 workstation and make veterinary dentistry a rewarding part of your day.

We have had the occasion to call on your service staff for advice and repairs regarding our dental machine. In every case, your staff have been very knowledgeable and provided assistance with patience and courtesy. You should be congratulated for outstanding service to your clients and we will continue to deal exclusively with you for our dental machine requirements."

Peter Roche,

Ferny Grove Veterinary Surgery

7. LEGAL INFORMATION

EMC - ELECTROMAGNETIC COMPATIBILITY

EMC is the interaction of electrical equipment with its electromagnetic environment and / with other equipment, ensuring the electronic devices operate with safe EMI levels and prevents interference that can potentially lead to malfunctioning, equipment damage, or even equipment failure.

Evolution Professional Veterinary Dental Unit:

EN IEC 55014-1 : 2021

EN IEC 55014-2 : 2021

EN IEC 61000-6-1 : 2019

EN IEC 61000-6-4 : 2019

FCC CFR 47 Part 15.109

EN IEC 61000-3-2: 2019 + A1 : 2021

EN 61000-3-3: 2013 + A1 : 2019

CE - EUROPEAN CONFORMITY

CE certification is a declaration of compliance with the European standards for health, safety and environmental protection for products sold or moved within the European Economic Area or other countries that recognize it.

UKCA - UK CONFORMITY ASSESSED

UKCA marking indicates conformity with the applicable requirements for products sold within Great Britain. The UKCA marking is used for products that meet EU requirements in Great Britain, Northern Ireland and the EU.

8. SPARE PART & SERVICE LIST

AS001- ADVANCED SERVICE KIT USED EVOLVE

- » **PS001** – 1 iM3 straw; Dental unit water treatment
- » **P7957** – 1 Exhaust bottle bacterial filter complete – Purple
- » **I64840** – 1 Oil free compressor yellow motor intake filter
- » **P1414** – Green tubing.
- » **W0530** (Tubing) – 1 Filter for water bottles and ultrasonic scalers.

AS002- ADVANCED SERVICE SUCTION KIT USED ON ADVANCED AND PROFESSIONAL

- » **PS001** – 1 iM3 straw; Dental unit water treatment
- » **P7957** – 2 Exhaust bottle bacterial filter complete – Purple.
- » **I64840** – 1 Oil free compressor yellow motor intake filter.
- » **P1414** – 1 Green tubing
- » **W0530** (Tubing) – Filter for water bottles and ultrasonic scalers
- » **MD30** – iM3 Advantage Handpiece Oil -30ml
- » **MD400** – Advantage Handpiece oil Aerosol
- » **PS001** – iM3 Straw
- » **I6580** – CLS / Water Toggle Switch

24 MONTHS WARRANTY

iM3 is liable for material or manufacturing defects within a warranty period of 24 months from the date of purchase.

We accept no responsibility for damage caused by incorrect handling or repairs carried out by the third parties not authorized to do so by iM3!

Claims under warranty- accompanied by proof of purchase- must be sent to the vendor or to an authorized iM3 service organization . The provision of service under warranty extends neither the warranty period nor any other guarantee period.

9. TROUBLE SHOOTING & MAINTENANCE

COMPRESSOR

 Isolate power supply before opening the compressor enclosure

IDENTIFYING LEAKS

To ensure the Compressor Motor is not subjected to higher than normal operating temperatures, each iM3 Dental unit should be checked for air leaks at regular intervals, the easiest way to check for air leaks is to pump the compressor tank up to 110 psi. / 7.5 BAR. and turn the motor off.

Make sure the airline to the dental head is connected with the pressure switch for the water bottles also turned on. Once the tank has been pumped to 110 psi. / 7.5 BAR., check the pressure in the tank 10 minutes later.

If the pressure has dropped more than 10% during that time then the unit is leaking air which should be identified and remedied to ensure the compressor motor does not overheat due to excessive running caused by the air leak.


Fig.46 - Inner workings of Dental Unit

ADJUSTING REGULATOR AIR PRESSURE

With the Dental Unit pressurised and air lines turned ON, the black knob on top of the regulator is pulled up from the locked to the free position to adjust air pressure (see picture page 21). Depress the black knob after setting to lock pressure adjustment. The regulator gauge registers working or air line pressure which may be increased by turning clockwise and lowered by turning anti clockwise. The pressure air line gauge will register the working pressure, which should be 70 psi./ 6BAR. Please note the regulator pressure us preset by uM3 before delivery of your unit

CIRCUIT BREAKER PUSH TO RESET

Unit starts automatically. Isolate power supply and relieve air pressure by opening tank drain valve prior to commencing non electrical service tasks

Danger High voltage used on this equipment. Isolate from mains before removing any covers, which should only be carried out by a qualified electrician . Ensure power outlet is adequately rated for this unit.

COMPRESSOR MAINTENANCE

FEATURE	DAILY	WEEKLY	MONTHLY	YEARLY
Compressor		Turn off power and unplug from wall. Drain moisture from main compressor tank by opening the tank drain on the top of the white tank marked with a yellow sticker. Catch the fluid from the white drain tube into a container. Keep the drain open until all air is removed from the tank.		
Regulator Auto Drain Filter	Automatic when tank drained of moisture			Replace regulator filter
Regulator Air Pressure	Check that correct air pressure is being delivered			

For a yearly service kit or to arrange for a yearly service of your dental unit please contact support@im3dental.com or www.im3vet.com

HANDPIECE MAINTENANCE

Lubrication of HS & LS Handpieces

Refer to separate handpiece manual

All handpieces should be lubricated regularly. We recommend the use of iM3 MD-30 handpiece oil before use and after 30 minutes of running time. The foot pedal should be depressed for 30 seconds after oiling to ensure the lubricant will reach the turbines and working portions of the handpiece.

Don't forget to depress the foot control after oiling to ensure lubricant reaches the turbines

All Handpieces have a working life which can be extended by maintaining the correct operating pressures and speeds as well as adequate lubrication.

The working life of most Handpieces turbines is 12-24 months, however, this will depend upon the degree of care taken and the amount of use.

Preventative maintenance will ensure a working Handpiece at all times.


Fig.47 - Lubrication of Handpieces

HANDPIECE	DAILY	WEEKLY	MONTHLY	YEARLY
iM3 Low Speed Handpiece Oil in smaller of 2 big holes. Depress the foot pedal for 10 seconds to ensure the oil reaches the Turbine	Oil AM or PM before or after use. 2-3 drops Use iM3 MD-30 Lubricant.			Replace turbine and 'o' rings and sealing gasket. If worn
iM3 High Speed Handpiece as per the instructions for the Low speed above.	Oil AM or PM before or after use. 2-3 drops Oil with MD-30 Lubricant.			Replace turbine and 'o' rings and sealing gasket. If worn
Prophy Head R&R Oscillating	Dispose of prophy head between patients to prevent cross contamination			

HANDPIECE TROUBLE SHOOTING

PROBLEM	CAUSE	SOLUTION
No Water Flow in Handpiece	1. Compressor not ON.	1. Turn compressor ON and check that ON/OFF tap is in the ON position.
	2. Water Flow Control not ON (P4- E).	2. Turn water control to ON.
	3. Water/CLS Bottle Empty.	3. Fill water/CLS Bottle(s).
	4. Check for a 'kink' in the water line.	4. Check lines and pressure.
	5. Check if 3 way syringe has Water.	5. If yes, lack of water must be in the valves or the handpiece. The 3 way syringe is the quickest way to check for water.
	6. Use of tap water, blocking water valves.	6. Valves will need replacing.
	7. Low regulator pressure	7. Check regulator is set at 70 PSI.
Mineral Accumulation	1. Use of tap water - mineral or chemical action blockage will usually take 3-4 years to occur.	1. Pass a fine wire through the water intake hole. Pass another fine wire through the hole in the handpiece head. Pass the wire as far as possible each way. 'Blow out' the loosened mineral deposits.
	2. Verdigris (Green Slime) build up in water valves due to chemical reaction with tap water.	2. Replace valves
Use of Chlorhexidine	1. Most brands of Chlorhex will deposit out of solution and block valves, etc. CLS is designed to be used in its diluted strength without depositing. NOTE: Do not dilute CLS beyond its label recommendations or it could result in deposits.	1. Contact your local iM3 Official Service Centre
	2. Chlorine in tap water will cause a flocculation and block filters and valves	2. Contact your local iM3 Official Service Centre. Use only distilled water

LOW SPEED TROUBLE SHOOTING

PROBLEM	CAUSE	SOLUTION
Handpiece won't run	1. No or low air pressure.	1. Check the air line to be sure the air is turned on and the line is unobstructed.
	2. Air line improperly sealed.	2. Reconnect air line to handpiece making sure the gasket is fully contracted.
	3. Inferior oil; too much or excess more commonly too little oil or water in the turbine.	3. See lubrication instructions on page 17.
	4. Broken drive head. Note new handpieces require a "Run in time" this can be up to 20 hours of use	4. Contact your local iM3 official Service Centre.
		If the motor stops spinning place a fine screw driver down the nose and rotate shaft.
Can't fit R&R prophyl angle	1. Plastic drive shaft broken off internally in Chuck housing.	1. Contact your local official iM3 Service Centre.
	2. Improper lubrication. Only fit plastic iM3 R&R prophyl angles Code L7466.	2. See lubrication instructions

CARTRIDGE TOGGLE

Toggle switches become floppy

All iM3 dental units use a unique cartridge system for the toggle switches. The toggles are used to pressurize the water/CLS bottles. When the o-rings become worn you simply push the cartridge in and turn a 1/4 turn anti clockwise. The cartridge can then be removed from the valve casing. To install a new toggle, locate the desired position, push and rotate a 1/4 turn clockwise. New toggles can be ordered from iM3 and installed by the user. Toggles can be installed in the 12, 3, 6 and 9 o'clock positions. If the toggle is in the wrong position simply remove and reposition. The Dental unit should be turned off and depressurised before replacement


Fig.48 - Toggle switch for Water and CLS

ULTRASONIC SCALER HANDPIECE

 For cleaning information read in full the sub manual for the Scaler

10. TECHNICAL PRODUCT DETAILS

iM3 is liable for material or manufacturing defects within a warranty period of 24 months from the date of purchase.

We accept no responsibility for damage caused by incorrect handling or repairs carried out by the third parties not authorized to do so by iM3!

Claims under warranty – accompanied by proof of purchase – must be sent to the vendor or to an authorized iM3 service organization. The provision of service under warranty extends neither the warranty period nor any other guarantee period.

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