



Low Lift Scissor Table II

# More comfort for your back!



604200

## Low Lift Scissor Table II

The EICKEMEYER® Low Lift Scissor Table II is perfect for lifting and transporting larger patients without unnecessary handling and movement. The table is available in 8 different colours and the optional add-ons make this table a real aa-rounder for your practice.

- ▶ Minimal height of 28 cm allowing large dogs to walk onto the table
- ▶ Maximum height: 105 cm
- ▶ Electric operating table helps protect your back
- ▶ Table height adjustable via foot pedal
- ▶ Stainless steel construction
- ▶ Smooth flat top with anti-drip edges and drain hole
- ▶ Wooden plate under table top to absorb vibration and reduce noise
- ▶ Maximum load: 100 kg
- ▶ Tabletop dimensions (in cm): L 130 x W 60
- ▶ Includes: Power Supply and Tie Down Cleats

Variations			
Colour	RAL value	Colour name	Item
		Stainless steel	<b>604200</b>
	3000	Flame Red	<b>604201</b>
	6018	Lime Green	<b>604202</b>
	5022	Night Blue	<b>604203</b>
	2001	Red Orange	<b>604204</b>
	6005	Moss Green	<b>604205</b>
	4010	Pink	<b>604206</b>
	9010	Pure White	<b>604207</b>

### EICKEMEYER®

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom  
T 020 8891 2007  
info@eickemeyer.co.uk

TELEPHONE 020 8891 2007

[www.eickemeyer.co.uk](http://www.eickemeyer.co.uk)



Low Lift Scissor Table II

veterinary technology for life  
**Eickemeyer®**



604303



604305



604301



604304

#### **Additional Options**

In addition to the stationary table, you can choose one or more of the following table options:

#### **Mobile Kit**

- ▶ Includes Wheels, Battery and Charger
- ▶ Table not included

**604306**

#### **Battery Kit**

- ▶ For wireless operation for stationary table
- ▶ Sold as battery only

**604307**

#### **Curved Dental Tray**

- ▶ Ergonomic, semi-circular removable drip tray with grid
- ▶ Made of stainless steel
- ▶ With removable tray with grid
- ▶ Tray size (in cm): W 50 x D 35

**604303**

#### **Integrated Scale Option**

- ▶ Weighing range: 0–120 kg
- ▶ ATTENTION: Cannot be combined with Tilt Mechanism (item no. 604301)

**604305**

#### **Electric Tilt Option**

- ▶ Foot pedal operated
- ▶ Tilt mechanism up to 35°
- ▶ Cannot be combined with Scale Option (item no. 604305)
- ▶ ATTENTION: Table not included

**604301**

#### **DIN Table Rail System**

- ▶ Lateral
- ▶ Attachable

**604304**

#### **EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom  
T 020 8891 2007  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**