



KOWA SL-19  
Portable Slit Lamp

veterinary technology for life  
**Eickemeyer**®

# Look through my eyes



## **KOWA SL-19 Portable Slit Lamp**

Now more ergonomic and even lighter! With its rounded shape and curved handle design, the new KOWA SL-19 Slit Lamp is even more ergonomic. Furthermore, the weight of the hand-held unit has been significantly reduced.

### **New features**

- ▶ New LED backlight that can be used simultaneously with the slit light
- ▶ New LED blue light illumination for fluorescence observation
- ▶ Light intensity can be adjusted with one hand
- ▶ Infinitely adjustable slit length (1.5 – 12 mm)
- ▶ Longer operating time: up to 6 hours at maximum intensity
- ▶ Available in 5 colors: white, black, blue, green and pink

### **Technical data**

- ▶ Magnification: 10x and 16x
- ▶ Working distance: about 80 mm at 16x magnification, about 100 mm at 10x magnification
- ▶ Effective field of view: Ø 10 mm at 16x magnification, Ø 15mm at 10x magnification
- ▶ Slit width: Ø 0.1, Ø 0.2, Ø 0.8 mm
- ▶ Spot size: Ø 1, Ø 5, Ø 12 mm (slit length continuously adjustable from 1.5 – 12 mm)
- ▶ Light intensity: 200 – 20,000 Lux (continuously adjustable)
- ▶ Power supply: 3x AA batteries (DC 4.8 – 6.4 V)
- ▶ Dimensions (in mm): W 197 x D 107 x H 238
- ▶ Weight (without batteries): only 620 g

**173690**

### **EICKEMEYER®**

3 Windmill Business Village  
Brooklands Close, Sunbury-on-Thames  
Surrey, TW16 7DY, United Kingdom

T 020 8891 2007  
info@eickemeyer.co.uk

**TELEPHONE 020 8891 2007**

**www.eickemeyer.co.uk**



KOWA SL-19  
Portable Slit Lamp

### Optional accessories for KOWA SL-19

- ▶ For digital image documentation using smartphones
- ▶ Affordable, efficient and powerful accessory for recording videos and images
- ▶ Connect or remove within seconds
- ▶ Smartphone not included

#### iPhone Adapters

iPhone 14 Pro Max Adapter  
**173635**

iPhone 12 Pro Max Adapter  
**176341**

iPhone 14 Pro Adapter  
**173634**

iPhone 12 Pro Adapter  
**176343**

iPhone 14 Plus Adapter  
**173633**

iPhone 12 Adapter  
**176340**

iPhone 14 Adapter  
**173632**

iPhone 12 mini Adapter  
**173642**

iPhone 13 Pro Max Adapter  
**173637**

iPhone 11 Pro Max RP Adapter  
**173666**

iPhone 13 Pro Adapter  
**173638**

iPhone 11 Pro RP Adapter  
**173665**

iPhone 13 Adapter  
**173639**

iPhone 11 RP Adapter  
**173664**

iPhone 13 mini Adapter  
**173636**

#### Samsung Adapters

Samsung Galaxy S10+ Adapter  
**173644**

Samsung Galaxy S9+ Adapter  
**173657**

Samsung Galaxy S10 Adapter  
**173658**

Samsung Galaxy S9 Adapter  
**173656**

Samsung Galaxy S10e Adapter  
**173659**

