

iM3 - The Veterinary Dental Company

Complete Tungsten Carbide Bur Range from iM3

Bur Type		Head Size (0.1 mm)	Cutting Length (mm)	Friction Grip for High Speed Handpieces			SW Handpieces
				FG Feline (16 mm)	FG Standard (19 mm)	FG Surgical (25 mm)	RA / Latch (22 mm)
Round / Rose	1/4	005	-	D7910	D7911	D7915	C9091
	1/2	006	-	D7909	D7912	D7916	C9092
	1	008	-		C9050		C9093
	2	010	-	D7908	C9080	C9085	C9095
	3	012	-				C9096
	4	014	-	D7920	C9110	C9115	C9100
	5	016	-				C9097
	6	018	-	D7907	C9140	C9145	C9160
Cross-cut Fissure Taper	8	023	-	D7930	C9142	D7935	C9094
	699	009	4.2		D8012		
	700	010	4.2	C9001	D8013	D8015	C9098
	701	012	4.2	C9002	C9320	D8016	C9330
	701L	012	6.0		C7010	C7015	
	702	016	4.4		C9350		C9390
	702L	016	6.0		D7810	D7815	
	703	021	4.8		C7030		
Cross-cut Fissure	703L	021	7.5		C4450		
	556	009	4.2		C9279	C9275	
	557	010	4.2		C9280		C9290
	557L	010	6.0		C9278		
	558	012	4.2		C9250	C9255	
Fissure Taper	559	014	4.4		C9300		
	169	009	4.2		D8011		
	169L	009	5.2		D8006	D8014	
	170L	010	6.0		D8009	C9089	
	171L	012	6.0		D8008		
Round End Xcut Fissure Taper	172L	016	6.0		D8007		
	1701	012	4.2		D7949		
Round End Xcut Fissure	1702	016	4.4		D7950		
	1557	010	4.2			D7801	
Round End Fiss.	1558	012	4.2			D7802	
	1158	012	4.2				
Inverted Cone	33 1/2	006	0.5	D7940	C9229	D7803	
	35	010	1.1		C9230		
	37	014	1.4		C9231		D7830
Pear	329	006	1.2		D8329		
	330	008	1.6		C9440		
	331	010	1.7		D8331		
	332L	012	4.2		D8332		
	333L	014	4.4		D8333		
Others	Arkansas Flame	025	7.0		D7851		D7850
	Metal Cutter 557E	010	4.2		D8557		
	Flame 7903	012	3.6		D8010		
	Lindemann 166	021	10.0				D8559
	Lindemann 162	016	9.0		D8558		



Round / Rose



Cross-cut Fissure Taper



Cross-cut Fissure Taper L



Cross-cut Fissure



Cross-cut Fissure L



Fissure Taper



Fissure Taper L



Flame / Fluted



Round End Xcut Fiss Taper



Round End Xcut Fissure



Round End Fissure



Inverted Cone



Pear



Metal Cutter



Arkansas Flame



Lindemann