

Transducers for **Vetus E5**



Probe	Illustration	Type of Array	Parameters	Application	Biopsy kit
C6-1m		Curved	- Bandwidth: 1.2-6.0 MHz - Extended FOV: 74° - Convex Radius: 55mm	Abdomen (Canine, Equine, Bovine, Ovine), Reproduction (Canine, Equine, Bovine, Ovine), MSK (Canine, Equine, Bovine, Ovine)	NGB-055, multi angle, reusable
mC11-3m		Curved	- Bandwidth: 2.6-12.8 MHz - Extended FOV: 114° - Convex Radius: 15mm	Abdomen (Canine, Feline), Cardiology (Canine, Feline), Reproduction (Canine, Feline), Small Parts (Canine, Feline)	NGB-018 multi angle, reusable
L14-3m		Linear	- Bandwidth: 3.0-14.0 MHz - FOV (max) - Non-extended: 3.8cm	Abdomen (Canine, Feline), Small Parts (Canine, Feline), MSK (Canine, Feline, Equine)	NGB-007, multi angle, reusable
L13-3m		Linear	- Bandwidth: 3.0-13.0 MHz - FOV (max) - Non-extended: 3.8cm	Abdomen (Canine, Feline), Small Parts (Canine, Feline), MSK (Canine, Feline, Equine)	NGB-007, multi angle, reusable

